

# LAARMANN<sup>®</sup>

*Innovators in Solids*

## PLASTICS CM1000



Laarmann Group B.V.  
Op het Schoor 6  
6041 AV Roermond  
Netherlands

[info@laarmann.eu](mailto:info@laarmann.eu)  
[www.laarmann.eu](http://www.laarmann.eu)  
[webshop.laarmann.eu](http://webshop.laarmann.eu)  
+31 6 23 40 00 33

# PLASTICS APPLICATION

## Grinding parameters

Grinding system: Cutting Mill Type CM1000

Accessories: Bottom sieve 5,00mm

Sample material: Plastics

Product feed size: Pieces of up to 80mm length

Grinding time: 90 seconds

Sample amount: 180 gram (see picture 3)

76,73 gram (see picture 4)

41,30 gram (see picture 5)



## Recommendation

For the sample preparation of a.m. plastic materials we recommend our Cutting Mill Type CM1000 under the above mentioned conditions.